# Instructions For Use: Osteon Scan Gauge Kit – GEN 1.0 – All Jurisdictions

ART-012, Version 5.0 (Effective from 21/12/2023)

en

## Indication for use:

Nexus iOS Scan Gauge Kits are designed to facilitate the acquisition of accurate implant position readings in scans taken with intraoral scanners for full arch prosthetic restorations. The scan gauges are to be mounted directly on the implant abutment embedded in the patient's mouth using the driver and screws provided in the kit.

## Validated Intraoral Scanners:

Below are the Intraoral Scanners that have been validated for use with the Nexus iOS Scan Gauge System:

Brand Name	Intraoral Scanner Name
3Shape	TRIOS 3
3Shape	TRIOS 4
Dentsply Sirona	Primescan
Dentsply Sirona	Cerec Primescan
iTero	Element
iTero	Element 2
iTero	Element 5D
MEDIT	i500
MEDIT	I700
Shining 3D	Aoralscan 3D
AlliedStar	AS 200E

## Contraindication:

Titanium (Grade 5) alloys used in Scan Gauges and screws are biocompatible; however, some patients may present allergies or hypersensitivity.

## Warnings and Precautions:

- Mounting of scan gauges onto the patient's mouth should only be conducted by dental professionals.
- Ensure the seating surface of the implants are clean and free of residue before securing the scan gauges.
- The scan gauges are small and should be handled with caution to avoid the patient from swallowing it.
- Each Scan gauge kit is compatible with one or few specific implant platforms. Using the scan gauge kit on a noncompatible implant platform would generate inaccurate scans and is not recommended.
- Scan gauges should always be secured with screws compatible with the respective implant system. Screw compatibility information can be accessed from the website.
- Any modifications (grinding, deburring, or otherwise retouching) to the scan gauge implant interface or the seven facets of the scan gauges would damage the scan gauge and result in inaccurate scans.
- 18 scan gauges provided within a kit is unique to the kit and should not be mixed or replaced with scan gauges provided with another kit.
- Refer to scanning instructions specific to each scanner provided with the kit to get a more thorough

understanding of the scanning process required for this system.

## Cleaning:

Cleaning and decontamination of the kit and its contents must be performed before and after every use to ensure removal of most contaminants resulting from medical activity. The following methods are recommended – thermal disinfection (90°C for 5 minutes) or semi-automatic ultrasonic cleaning (15 minutes) with a washer-disinfector of neutral pH. The box could be rinsed under flowing water while brushing with a soft brush to remove any residue. The brush should not damage the surfaces of the gauge as it would make them non-functional. Following cleaning, the kit and its contents could be dried using clean compressed air.

## Sterilization:

The scan gauge kit and its contents are supplied non-sterile and MUST be sterilized prior to use. Remove the foam insert prior to using the kit. Before sterilization remove screw from packaging and place inside the screw tray provided within the kit

The following sterilisation options are recommended:

1. Cycle type: Steam sterilization Pre-Vacuum Cycle, exposure time: 4 minutes, temperature: 134°C (273°F), drying time: 40 minutes.
2. Cycle type: Steam sterilization Gravity Cycle, exposure time: 20 minutes, temperature: 134°C (273°F), drying time: 40 minutes

Also, ensure there is no contact or overlap between the contents in the gauge as poor positioning may result in insufficient sterilisation.

## Safety information:

Osteon products have not been evaluated for safety when exposed to electric, radiation and thermal energy and compatibility is unknown. Exposure to these energies may result in patient injury and adequate precautions should be taken.

**Note:** Any serious incident(s) that has occurred in relation to the device shall be reported to Osteon and the competent authority having jurisdiction in the locale where the incident occurred.

## MRI Safety Information:

Warning: The RF safety of the device has not been tested.

## Storage:

Ensure the kits stored in a way that protects it from excessive temperatures and humidity.

Product of Australia

Worldwide: +61 3 9264 0111 | Australia: 03 9264 0111 | [info@osteonmedical.com](mailto:info@osteonmedical.com) | [www.osteonmedical.com](http://www.osteonmedical.com)  
759-767 Springvale Rd, Mulgrave VIC 3170 Australia

# Instructions For Use: Osteon Scan Gauge Kit – GEN 1.0 – All Jurisdictions

ART-012, Version 5.0 (Effective from 21/12/2023)

## Labelling Symbols:

Labelling Symbols according to ISO 15223-1:2021	
Symbol	Description
	Manufacturer
	Date of manufacture
	Non-Sterile
	Consult electronic Instructions for Use (E-IFU)
	Fragile, handle with care
	Medical Device
	CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
	MR Unsafe
	Authorized representative in the European Union
	Catalogue number
	Serial Number
	Importer
	Distributor
	Do not use if package is damaged
	Keep dry (before opening)

	Unique Device Identifier
	Responsible Person (United Kingdom)

## Disposal:

For disposing product and packaging comply with applicable national waste disposal regulations in your country.

## Contact Details:

**Australia:**  
**Implant Solutions Pty Ltd (trading as Osteon Medical)**  
 759-767 Springvale Rd, Mulgrave, VIC 3170, Australia  
**Contact:**  
[customerservice.au@osteonmedical.com](mailto:customerservice.au@osteonmedical.com)  
 +61 3 9264 0111

**USA:**  
**Osteon Medical North America**  
 Suite 150, 175 Technology Drive, Irvine, CA  
**Contact:**  
[customerservice.USA@osteonmedical.com](mailto:customerservice.USA@osteonmedical.com)  
 1-888-203-6180

**Europe :**  
**Osteon Medical Europe**  
 11A rue des Aulnes, Champagne au Mt D'or, 69410  
**Contact:**  
[customerservice.eu@osteonmedical.com](mailto:customerservice.eu@osteonmedical.com)  
 06 01 46 33 37

**Canada:**  
**Keystone Dental Group**  
 154 Middlesex Turnpike, Burlington MA 01803, USA  
**Contact:**  
[KPranjivan@keystonedental.com](mailto:KPranjivan@keystonedental.com)  
 (416) 660-1727

## Authorised Representatives:

**EC REP** **European Healthcare & Device Solutions (Ireland) Ltd.**  
 Stratton House, Bishopstown Road, Bishopstown, Cork T12 Y9TC, Ireland  
[tqmuk@aol.com](mailto:tqmuk@aol.com)

**UK RP** **European Device Solutions Ltd.**  
 15 Coanwood Drive, Whitley Bay, Tyne & Wear NE25 9GB, UK  
[Prowden2510@gmail.com](mailto:Prowden2510@gmail.com)

**US Agent** **Keystone Dental Group**  
 154 Middlesex Turnpike, Burlington, MA US 01803  
[NDeAngelo@keystonedental.com](mailto:NDeAngelo@keystonedental.com)

Product of Australia

Worldwide: +61 3 9264 0111 | Australia: 03 9264 0111 | [info@osteonmedical.com](mailto:info@osteonmedical.com) | [www.osteonmedical.com](http://www.osteonmedical.com)  
759-767 Springvale Rd, Mulgrave VIC 3170 Australia

# Instructions For Use: Osteon Scan Gauge Kit – GEN 1.0 – All Jurisdictions

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