en

Indication for use:

Nexus iOS Scan Gauge Kits are designed to facilitate the acquisition of accurate implant position readings in scans taken with intraoral scanners for full arch prosthetic restorations. The scan gauges are to be mounted directly on the implant abutment embedded in the patient's mouth using the driver and screws provided in the kit.

Contraindication:

Titanium (Grade 5) alloys used in Scan Gauges and screws are biocompatible; however, some patients may present allergies or hypersensitivity.

Warnings and Precautions:

- Mounting of scan gauges onto the patient's mouth should only be conducted by dental professionals.
- Ensure the seating surface of the implants are clean and free of residue before securing the scan gauges.
- The scan gauges are small and should be handled with caution to avoid the patient from swallowing it.
- Each Scan gauge kit is compatible with one or few specific implant platforms. Using the scan gauge kit on a noncompatible implant platform would generate inaccurate scans and is not recommended.
- Scan gauges should always be secured with screws compatible with the respective implant system. Screw compatibility information can be accessed from the website.
- Any modifications (grinding, deburring, or otherwise retouching) to the scan gauge implant interface or the seven facets of the scan gauges would damage the scan gauge and result in inaccurate scans.
- 18 scan gauges provided within a kit is unique to the kit and should not be mixed or replaced with scan gauges provided with another kit.
- Refer to scanning instructions specific to each scanner provided with the kit to get a more thorough understanding of the scanning process required for this system.

Cleaning:

Cleaning and decontamination of the kit and its contents must be performed before and after every use to ensure removal of most contaminants resulting from medical activity. The following methods are recommended – thermal disinfection (90°C for 5 minutes) or semi-automatic ultrasonic cleaning (15 minutes) with a washer-disinfector of neutral pH. The box could be rinsed under flowing water while brushing with a soft brush to remove any residue. The brush should not damage the surfaces of the gauge as it would make them non-functional. Following cleaning, the kit and its contents could be dried using clean compressed air.

Sterilization:

The scan gauge kit and its contents are supplied non-sterile and MUST be sterilized prior to use. Remove the foam insert prior to using the kit. Before sterilization remove screw from packaging and place inside the screw tray provided within the kit. Any of the following sterilization method could be

employed - Steam sterilization (in accordance with ISO 17665 or other national standard), Moist Heat sterilization, Gamma Radiation, Ethylene sterilization or Hydrogen peroxide sterilization. The following is recommended – cycle type: saturated steam, exposure time: 18 minutes, temperature: 134°C (273°F), drying time: minimum 20 minutes. Also, ensure there is no contact or overlap between the contents in the gauge as poor positioning may result in insufficient sterilisation. Flash sterilization should not be performed. The manufacturer disclaims any liability for sterilization performed by the end user that is non-compliant with the ISO 17665 standard.

Safety information:

Osteon products have not been evaluated for safety when exposed to electric, radiation and thermal energy and compatibility is unknown. Exposure to these energies may result in patient injury and adequate precautions should be taken.

Note: Any serious incident(s) that has occurred in relation to the device shall be reported to Osteon and the competent authority having jurisdiction in the locale where the incident occurred.

MRI Safety Information:

Warning: The RF safety of the device has not been tested.

Storage:

Ensure the kits stored in a way that protects it from excessive temperatures and humidity.

Labelling Symbols:

Labelling Symbols according to ISO 15223-1:2021	
Symbol	Description
	Manufacturer
	Date of manufacture
	Non-Sterile
	Consult electronic Instructions for Use (E-IFU)
	Fragile, handle with care
	Medical Device

Instructions For Use: Osteon Scan Gauge Kit – GEN 1.0 – All Jurisdictions

ART-012, Version 3.0 (Effective from 11/01/2023)

	CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
	MR Unsafe
	Authorized representative in the European Union
	Catalogue number
	Serial Number
	Importer
	Distributor
	Do not use if package is damaged
	Keep dry (before opening)
	Unique Device Identifier
	Responsible Person (United Kingdom)

Contact Details:

Australia:
Implant Solutions Pty Ltd (trading as Osteon Medical)
 759-767 Springvale Rd, Mulgrave, VIC 3170, Australia
Contact:
customerservice.au@osteonmedical.com
 +61 3 9264 0111

USA:
Osteon Medical North America
 Suite 150, 175 Technology Drive, Irvine, CA
Contact:
customerservice.USA@osteonmedical.com
 1-888-203-6180

Europe :
Osteon Medical Europe
 11A rue des Aulnes, Champagne au Mt D'or, 69410
Contact:
customerservice.eu@osteonmedical.com
 06 01 46 33 37

Canada:
Keystone Dental Group
 154 Middlesex Turnpike, Burlington MA 01803, USA
Contact:
KPranjivan@keystonedental.com
 (416) 660-1727

Authorised Representatives:

EC REP **European Healthcare & Device Solutions (Ireland) Ltd.**
 Stratton House, Bishopstown Road, Bishopstown, Cork T12 Y9TC, Ireland
tgmuk@aol.com

UK RP **European Device Solutions Ltd.**
 15 Coanwood Drive, Whitley Bay, Tyne & Wear NE25 9GB, UK
Prowden2510@gmail.com

US Agent **Keystone Dental Group**
 154 Middlesex Turnpike, Burlington, MA US 01803
NDeAngelo@keystonedental.com

Disposal:

For disposing product and packaging comply with applicable national waste disposal regulations in your country.

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