# Instructions For Use: Osteon Precision Milled CoCr Suprastructure – AU, EU, and UK

ART-020, Version 2.0 (Effective from 03/11/2022)

en

## Indication for use:

The Osteon Precision Milled Suprastructure is indicated for attachment to dental abutments in the treatment of partially or fully edentulous jaws for the purpose of restoring chewing function.

## Description:

Each suprastructure, made for an individual patient as per the specifications provided by the prescribing health professional, is digitally designed using advanced computer aided design (CAD) system, precision milled (CAM) in biocompatible materials is compliant with the requirements of relevant Essential Principles. The device is to be retained on implants with screws.

## Material:

The bar is made of Cobalt Chrome Alloys as per DIN EN ISO 22674, free of Beryllium and Nickel. The overlay is fabricated of Acrylic (PMMA) or Ceramic (Zirconia).

## Contraindications:

Although material used is biocompatible, some patient may present allergies or hypersensitivity to the alloy and its constituents, which may include any or all of the following: Cobalt, Chrome, Tungsten, Silicon, Carbon, Iron, Manganese, Zirconia (Zirconia, Yttrium Oxide, Aluminium Oxide and other oxides), PMMA.

Patients with bruxism since an overload of the device may occur.

Cylinders are not to incorporate divergence that will cause the resulting angulation of any single construct to exceed the maximum angulation of the compatible multi-unit abutment system.

## Warnings and Precautions:

Proper surgical procedures and restorative techniques are the responsibility of the dental professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. Clinicians are also responsible for informing patients on proper maintenance of the device post restoration after evaluating their underlying health conditions and oral parafunction.

Small components dropped in the patient’s mouth could be swallowed or aspirated and therefore should be handled with care.

If any modifications are made to the implant suprastructure interface, the suprastructure may have a fitting issue. These products are made specifically for a patient. Re-use of the product for another patient should not be attempted. Osteon assumes no responsibility for attempted re-use of the device on another patient or defects that arise due to the modification (grinding, deburring, or otherwise retouching) of the device.

Potential causes of suprastructure fracture include but are not limited to inadequate implant support when attached to periodontally compromised teeth, non-passive fit of suprastructure, overloading due to improper occlusion,

incomplete seating of cemented abutments, and excessive cantilevering of pontics.

Divergence/Angulation: Cylinders are not to incorporate additional divergence that will cause the resulting angulation of any single construct to exceed the maximum angulation of the compatible multi-unit abutment system as stated in the Indications for Use of the Osteon Precision Milled Suprastructure.

## Sterilization:

Osteon products are supplied NON-STERILE and must be sterilized and/or disinfected according to a validated method prior to use in the oral environment. The recommended sterilization cycle for metallic bars is standard prevacuum steam sterilization, exposure at 132°C (270°F) for 4 minutes (after sterilizer has reached indicated temperature) with a dry time of 45 minutes, using a sterilization pouch that is cleared for the indicated cycle. For the overlays or a fully finished and cemented product (composed of bar and overlay), the clinician should use appropriate disinfection method prior to use in oral environment.

## Specific Instruction For Use:

Refer to the biaxial card provided for suprastructures containing biaxial systems to identify their precise location. Biaxial driver and OEM (original equipment manufacturer) driver should be used with biaxial and straight screws respectively. The Biaxial screws should be secured with the torque specified on biaxial screw packaging. For straight systems, screws should be torqued as per implant/abutment OEM specifications.

## Safety information:

Osteon products have not been evaluated for safety when exposed to electric, radiation and thermal energy and compatibility is unknown. Exposure to these energies may result in patient injury and adequate precautions should be taken.

**Note:** Any serious incident(s) that has occurred in relation to the device shall be reported to Osteon and the competent authority having jurisdiction in the locale where the incident occurred.

## MRI Safety Information:

Warning: The Osteon Precision Milled Suprastructures has not been evaluated for safety and compatibility in the MR Environment. It has not been tested for heating, migration, or image artefacts in the MR Environment. The safety of suprastructures in the MR Environment is unknown. Scanning a patient, who has this device may result in patient injury.

## Labelling Symbols:

Labelling Symbols according to ISO 15223-1:2021	
Symbol	Description
	Manufacturer
	Date of manufacture

Product of Australia

# Instructions For Use: Osteon Precision Milled CoCr Suprastructure – AU, EU, and UK

ART-020, Version 2.0 (Effective from 03/11/2022)

	Non-Sterile
	Consult electronic Instructions for Use (E-IFU)
	Fragile, handle with care
	Medical Device
	CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
	MR Unsafe
	Authorized representative in the European Union
	Catalogue number
	Serial Number
	Importer
	Distributor
	Do not use if package is damaged
	Keep dry (before opening)
	Single patient (multiple use)
	Unique Device Identifier
	Responsible Person (United Kingdom)

## Contact Details:

### Implant Solutions Pty Ltd (trading as Osteon Medical)

759-767 Springvale Rd, Mulgrave, VIC 3170, Australia

**Contact:**

[customerservice.au@osteonmedical.com](mailto:customerservice.au@osteonmedical.com)

+61 3 9264 0111

### Osteon Medical Europe

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**Contact:**

[customerservice.eu@osteonmedical.com](mailto:customerservice.eu@osteonmedical.com)

06 01 46 33 37

## Authorised Representatives:

### European Healthcare & Device Solutions (Ireland) Ltd.

Stratton House, Bishopstown Road, Bishopstown, Cork T12 Y9TC, Ireland

[tqmuk@aol.com](mailto:tqmuk@aol.com)

### European Device Solutions Ltd.

15 Coanwood Drive, Whitley Bay, Tyne & Wear NE25 9GB, UK

[Prowden2510@gmail.com](mailto:Prowden2510@gmail.com)

## Disposal:

For disposing product and packaging comply with applicable national waste disposal regulations in your country.

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06 01 46 33 37

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### EC REP

#### European Healthcare & Device Solutions (Ireland) Ltd.

Stratton House, Bishopstown Road, Bishopstown, Cork T12 Y9TC, Ireland  
[tqmuk@aol.com](mailto:tqmuk@aol.com)

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## Kontraindikation:

Although material used is biocompatible, some patient may present allergies or hypersensitivity to the alloy and its constituents, which may include any or all of the following: Cobalt, Chrome, Tungsten, Silicon, Carbon, Iron, Manganese, Zirconia (Zirconium Oxide, Aluminium Oxide and other oxides), PMMA.

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